

SPRINT[®]
PNS SYSTEM

Chronic Pain Relief



Life. Get
back to it.™


The breakthrough treatment
designed to provide sustained
relief from chronic pain.

Short-term treatment designed for long-term relief.

The SPRINT® System helps you get back to the activities you've missed most – without the worry of addiction or nerve damage.

60-day treatment designed to provide sustained pain relief:

- minimally-invasive, non-surgical
- drug-free pain treatment
- no permanent implant





After ten years of compartmentalizing the pain through meditation and yoga to just get through it, Patricia is now able to once again have a more rigorous, intense yoga practice. I'm "back to doing everything!...It's a miracle — the pain is gone."

— Patricia (patient with severe nerve pain related to Meralgia Paresthetica)

► Visit us at **SprintPNS.com**
to hear patient stories.



The SPRINT System

- The Pulse Generator sends gentle pulses through a very small, hair-like MicroLead™ to stimulate the nerve.
- The simple Hand-Held Remote allows you to easily control the level of stimulation.
- The majority of patients experience sustained relief after the MicroLead is removed and the treatment ends.





PULSE
GENERATOR



HAND-HELD
REMOTE

Clinically proven.

The SPRINT Peripheral Nerve Stimulation (PNS) System has been studied extensively for low back pain, shoulder pain, post-amputation pain, and post-operative pain, and is cleared for use up to 60 days.*

This breakthrough treatment offers **significant and sustained pain relief** and is the only stimulation device cleared that does not require incisions or a permanent implant.

MICROLEAD™
(ACTUAL SIZE)



Ask your doctor if
the SPRINT System
is right for you.

FAQ:

What can I expect?

You will have a minor procedure to implant the MicroLead(s). The SPRINT System will be left in place for up to 60 days, after which, the MicroLead(s) will be removed by your doctor.

Are there any activity limitations?

Expect to go about your normal activities while using the SPRINT System. Before taking a shower, you will need to disconnect the pulse generator. Bathing and swimming should be avoided during the treatment.

What risks are associated with SPRINT?

The most common risks are skin irritation and redness that usually go away within a few days without the need for medical care. The implanted MicroLead is constructed of very fine wire about the size of a human hair. It is possible that a lead remnant may remain beneath your skin following the procedure. No serious adverse events related to the use of the SPRINT System have been reported. For more information visit www.SprintPNS.com.

Insurance coverage

Many insurance companies cover the SPRINT System. SPRcare™ can provide support to you or your doctor's staff to obtain coverage from your insurance company.

Learn more at SprintPNS.com

SPRINT[®]
PNS SYSTEM

Breakthrough
Treatment for Pain

*The SPRINT PNS System is indicated for up to 60 days in the back and/or extremities for: (i) Symptomatic relief of chronic, intractable pain, post-surgical and post-traumatic acute pain; (ii) Symptomatic relief of post-traumatic pain; and (iii) Symptomatic relief of post-operative pain. The SPRINT PNS System is not intended to treat pain in the craniofacial region.

Physicians should use their best judgment when deciding when to use the SPRINT PNS System. For more information see the SPRINT PNS System IFU.

Most common adverse events are skin irritation and erythema. Results may vary. Rx only. View References for support of all claims at SprintPNS.com.

For information about Microlead Safety, visit us at <https://www.sprtherapeutics.com/patients/microlead-safety/>

The SPRINT[®] PNS System, Life. Get Back to It.[™], MicroLead[™], SPRcare[™] and SPR[®] are registered or common law trademarks of SPR Therapeutics, Inc. © SPR THERAPEUTICS, INC 2020.
9600-MKT-015-C

