

Prospective, Multicenter Study of Percutaneous Medial Branch PNS for the Treatment of Chronic Axial Low Back Pain

Christopher A. Gilmore, MD,¹ Mehul J. Desai, MD, MPH,² Thomas J. Hopkins, MD, MBA,³ Sean Li, MD,⁴ Michael J. DePalma, MD,⁵ Timothy R. Deer, MD,⁶ Steven P. Cohen, MD,⁷ Meredith J. McGee, PhD,⁸ Joseph W. Boggs, PhD⁸

¹Center for Clinical Research, ²International Spine, Pain & Performance Center, ³Duke University, ⁴Premier Pain Centers, ⁵Virginia iSpine Physicians, ⁶The Spine and Nerve Center of The Virginias, ⁷Johns Hopkins University, ⁸SPR Therapeutics, Inc.



Introduction

Chronic low back pain (LBP) is one of the most prevalent and challenging musculoskeletal conditions¹ and is the leading cause of disability in adults.

Minimally Invasive, 60-day Percutaneous PNS:

A promising non-opioid, non-destructive, and non-surgical treatment for chronic axial LBP, designed for use earlier in the treatment continuum than conventional neurostimulation.

Goal: Characterize responses to medial branch PNS in a prospective multicenter case series study in patients recalcitrant to multiple non-surgical treatments



Materials & Methods

Ongoing IRB-approved study; informed consent obtained from each subject.

Key Eligibility Criteria:



- Subjects with chronic axial LBP (≥ 3 months); no radicular pain
- Stable medication usage for ≥ 1 month prior to baseline
- No prior lumbar surgery or RFA within prior 6 months
- No anesthetic injections within prior 3 months
- Score of ≤ 20 on Beck Depression Inventory

PNS Lead Implantation: Bilateral, percutaneous PNS leads, targeting medial branches of the dorsal ramus in the center of pain

- Image Guidance:** ultrasound and/or fluoroscopy
- Confirmation:** Stimulation of medial branch confirmed by selective activation of multifidi

PNS Treatment: Stimulation for 6-12 hrs/day for up to 60 days

- Subjects continued normal activities
- Leads removed with gentle traction
- Long-term follow-up visits, up to 12 months after the 2-month PNS treatment (ongoing)



Results


Participant Demographics (n=74)

| | |
|--|-------------|
| Age (years) | 56.3 (13.5) |
| Body Mass Index (BMI) | 29.4 (4.6) |
| LBP Duration (years) | 16.0 (13.0) |
| Sex (% Female) | 53% |
| Baseline Opioid Usage (MME; n=20 on opioids at baseline) | 32.0 (37.1) |

Previously Failed LBP Treatments:


| | |
|----------------------------------|-----|
| Non-opioid Analgesics | 97% |
| Physical Therapy | 89% |
| Opioid Analgesics | 67% |
| TENS | 65% |
| Anesthetic or Steroid Injections | 57% |
| Epidural Injections | 46% |
| Radiofrequency Ablation | 23% |

LBP Diagnoses / Etiologies of Pain:



Percutaneous PNS Implantation


- 91% of participants received bilateral PNS leads
- L4 & L5 were the most commonly targeted spinal levels




Adverse Events: No serious or unanticipated adverse events (AEs).

- The most common AEs were mild skin irritation or pruritis (itching).


Reductions in Average Pain Intensity (BPI-5):



Reductions in Oswestry Disability Index (ODI):



Reductions in Pain Interference (BPI-9):




Majority of participants reported clinically² and statistically significant reductions in pain, disability, and pain interference with percutaneous PNS.

- Reductions were sustained through at least 14 months among a majority of participants.

Majority of participants reported statistically significant improvements in health-related quality of life (QoL)

- 91% reported QoL improvements with PNS:




- Statistically significant improvements in QoL with PNS (measured via RAND-36):

| Selected RAND-36 Subscales | Improvement |
|------------------------------------|---------------------|
| Role Limitations – Physical Health | 128%, 28.1 points** |
| Pain | 62%, 22.1 points** |
| Physical Functioning | 33%, 14.0 points** |
| Energy / Fatigue | 32%, 13.2 points** |

** p < 0.001

Reductions in Opioid Analgesic Usage:



Conclusions

- Clinically meaningful and statistically significant reductions in pain, disability, and pain interference were reported by a majority of participants who completed the Primary Endpoint and each of the long-term follow up visits.
- Given the minimally invasive, non-destructive nature of percutaneous PNS and the significant benefits, percutaneous PNS may provide a promising first-line neurostimulation treatment for patients with chronic axial LBP.

References:

¹ US Burden of Disease Collaborators. The State of US Health, 1990-2010 Burden of Diseases, Injuries, and Risk Factors. *JAMA*. 2013; 310(6):591-606.
² Dworkin et al. Interpreting the clinical importance of treatment outcomes in chronic pain clinical trials: IMMPACT recommendations. *J Pain* 2008; 9: 105-121.

Funding:

This study was funded by SPR Therapeutics, Inc.

